

**Precision Power International, Inc.**

25108 Marguerite Parkway

Suite A-159

Mission Viejo, CA 92692 USA

T 949 951 6784

F 949 916 6733

sales@precisionpowerinternational.com

Specification Sheet **LPS-AB3730XL-LT**
Version 20090508

Certified - Instrumentation Grade: 208/120 Vac 60 Hz cUL
Instrument Power Protection System (IPPS) for an Applied Biosystems PRISM 3730 or 3730xl Genetic Analyzer

Catalog Number: **LPS-AB3730XL-LT**

(1) Catalog Number LPS-AB3730XL-LT supports the following AB PRISM® 3700 series (208/120 Vac cUL) instrument elements:


	Desktop Workstation (Tower Computer & LCD Monitor)
AB PRISM® 3730	(1) 48 Capillary Genetic Analyzer or,
AB PRISM® 3730xl	(1) 96 Capillary Genetic Analyzer

PPI - Instrumentation Power Protection System (IPPS)
Specifications: ([complete system protection](#))

IPPS set with LPS & PDU combination



(IPPS shown relative size)



Applied Biosystems (AB) PRISM® 3730xl Genetic Analyzer (computer workstation not shown)

Instrument photo courtesy of Applied Biosystems

Design Installations: North America and Countries Utilizing NEMA Type Connectors (208/120 Vac 50/60 Hz)

Harmonization Code: 8504.90.4000

Export License: U.S. EAR (ECCN) EAR 99 - No License Required - NLR

Certification: IEC/IEEE 62040 Laboratory and Instrumentation Grade
UL 1778, CSA 22.2-107, IEEE 62.41 CAT. A Lightning Safety

Performance: Regenerated Power Technology (RPT), Digital Signal Processing (DSP)
True Sinewave Output, True Online, PWM Double Conversion (IGBT)
> 20,000 Hours MTBF